


Cannulated Screws / Staples

Lochschrauben / Staples

Tornillos canulados / Grapas


Viti cannulate / Staples

Vis perforée / Staples




**D-127-03 Cannulated Screws 3.5 mm and Instrument Set
196-045-17 Aluminium Case
196-127-05 Upper Tray
196-127-06 Lower Tray**

			Pieces
127-035-10 To 127-035-24	Cannulated Screw 3.5 mm (half thread) 10-24 mm	2 each.	16
127-035-26 To 127-035-40	Cannulated Screw 3.5 mm (half thread) 26-40 mm	4 each.	32
127-035-42 To 127-035-50	Cannulated Screw 3.5 mm (half thread) 42-50 mm	2 each.	10
127-045-10 To 127-045-24	Cannulated Screw 3.5 mm (full thread) 10-24 mm	2 each.	16
127-045-26 To 127-045-40	Cannulated Screw 3.5 mm (full thread) 26-40 mm	4 each.	32
127-045-42 To 127-045-50	Cannulated Screw 3.5 mm (full thread) 42-50 mm	2 each.	10
127-085-13	Washer 3.5 mm		6
127-160-27	Cannulated Drill Bit 2.7 mm Ø		1
127-160-35	Cannulated Drill Bit 3.5 mm Ø		1
127-168-35	Cannulated Tap 3.5 mm Ø		1
127-125-10	Guide Wire Ø 1.0 mm		1
127-175-35	Cannulated screwdriver		1
100-072-02	Screw Holding Forceps		10
127-212-35	Tap Sleeve 2.7 mm, drill sleeve 1.0 mm		1
127-200-35	Tap Sleeve 3.5 mm, drill sleeve 2.5 mm		1




D-127-04 Cannulated Screws 4.5 mm and Instrument Set
196-045-17 Aluminium Case
196-127-08 Upper Tray
196-127-09 Lower Tray

			Pieces
127-055-20 To 127-055-34	Cannulated Screw 4.5 mm (half thread) 20-34 mm	2 each.	16
127-055-36 To 127-055-44	Cannulated Screw 4.5 mm (half thread) 36-44 mm	4 each.	32
127-055-46 To 127-055-60	Cannulated Screw 4.5 mm (half thread) 46-60 mm	2 each.	16
127-065-20 To 127-065-34	Cannulated Screw 4.5 mm (full thread) 20-34 mm	2 each.	16
127-065-36 To 127-065-44	Cannulated Screw 4.5 mm (full thread) 36-44 mm	4 each.	20
127-065-46 To 127-065-60	Cannulated Screw 4.5 mm (full thread) 46-60 mm	2 each.	18
127-085-16	Washer 4.5 mm		6
127-160-32	Cannulated Drill Bit 3.2 mm \varnothing		1
127-160-45	Cannulated Drill Bit 4.5 mm \varnothing		1
127-168-45	Cannulated Tap 4.5 mm \varnothing		1
127-125-16	Guide Wire \varnothing 1.6 mm		1
127-175-45	Cannulated screwdriver		10
100-072-02	Screw Holding Forceps		1
100-115-45	Tap Sleeve for 4.5 mm Screw		1



D-127-07 Cannulated Screws 7.0 mm and Instrument Set
196-045-17 Aluminium Case
196-127-01 Upper Tray
196-127-02 Lower Tray

			Pieces
127-070-30 To 127-070-65	Cannulated Screw 7.0 mm \varnothing , 16 mm 30-35-40-45-50-55-60-65 mm	2 each.	16
127-070-70 To 127-070-92	Cannulated Screw 7.0 mm \varnothing , 16 mm 70-75-80-85-90-95-100 mm	4 each.	28
127-070-93 To 127-070-96	Cannulated Screw 7.0 mm \varnothing , 16 mm 105-110-115-120 mm	2 each.	8
127-080-45 To 127-080-65	Cannulated Screw 7.0 mm \varnothing , 32 mm 45-50-55-60-65 mm	1 each.	5
127-080-70 To 127-080-92	Cannulated Screw 7.0 mm \varnothing , 32 mm 70-75-80-85-90-95-100 mm	2 each.	14
127-080-93 To 127-080-96	Cannulated Screw 7.0 mm \varnothing , 32 mm 105-110-115-120 mm	1 each.	4
127-085-19	Washer 7.0 mm \varnothing		6
127-160-50	Drill Bit 5.0 mm \varnothing		1
127-168-70	Tap 7.0 mm \varnothing		1
127-175-70	Screw Driver		1
127-227-70	Tissue Protector Sleeve		1
127-125-20	Guide Wire 2.0 mm \varnothing 230 mm		10
100-072-20	Screw Holding Forceps		1



3.5 mm \varnothing with over flats S.W. 2.5 mm half / full threaded drill bit 2.7 mm

\varnothing 3.5mm, 4.5mm, 7.0mm



4.5 mm \varnothing with over flats S.W. 3.5 mm half / full threaded drill bit 3.2 mm


Dia of thread	3.5 mm	3.5 mm
Width over flats	2.5 mm	2.5 mm
Drill Bit	2.7 mm	2.7 mm
Length		
10 mm	127-035-10	127-045-20
12 mm	127-035-12	127-045-22
14 mm	127-035-14	127-045-24
16 mm	127-035-16	127-045-26
18 mm	127-035-18	127-045-28
20 mm	127-035-20	127-045-30
22 mm	127-035-22	127-045-32
24 mm	127-035-24	127-045-34
26 mm	127-035-26	127-045-36
28 mm	127-035-28	127-045-38
30 mm	127-035-30	127-045-40
32 mm	127-035-32	127-045-42
34 mm	127-035-34	127-045-44
36 mm	127-035-36	127-045-46
38 mm	127-035-38	127-045-48
40 mm	127-035-40	127-045-50
42 mm	127-035-42	127-045-52
44 mm	127-035-44	127-045-54
46 mm	127-035-46	127-045-56
48 mm	127-035-48	127-045-58
50 mm	127-035-50	127-045-60

Dia of thread	4.5 mm	4.5 mm
Width over flats	3.5 mm	3.5 mm
Drill Bit	3.2 mm	3.2 mm
Length		
20 mm	127-055-20	127-065-20
22 mm	127-055-12	127-065-22
24 mm	127-055-14	127-065-24
26 mm	127-055-16	127-065-26
28 mm	127-055-18	127-065-28
30 mm	127-055-20	127-065-30
32 mm	127-055-22	127-065-32
34 mm	127-055-24	127-065-34
36 mm	127-055-26	127-065-36
38 mm	127-055-28	127-065-38
40 mm	127-055-30	127-065-40
42 mm	127-055-32	127-065-42
44 mm	127-055-34	127-065-44
46 mm	127-055-36	127-065-46
48 mm	127-055-38	127-065-48
50 mm	127-055-40	127-065-50
52 mm	127-055-42	127-065-52
54 mm	127-055-44	127-065-54
56 mm	127-055-46	127-065-56
58 mm	127-055-48	127-065-58
60 mm	127-055-50	127-065-60



7 mm \varnothing with over flats SW 5.0 mm 16 mm length of thread, drill bit 5.0 mm

127-070-25	25 mm	127-070-75	75 mm
127-070-30	30 mm	127-070-80	80 mm
127-070-35	35 mm	127-070-85	85 mm
127-070-40	40 mm	127-070-90	90 mm
127-070-45	45 mm	127-070-91	95 mm
127-070-50	50 mm	127-070-92	100 mm
127-070-55	55 mm	127-070-93	105 mm
127-070-60	60 mm	127-070-94	110 mm
127-070-65	65 mm	127-070-95	115 mm
127-070-70	70 mm	127-070-96	120 mm



7 mm \varnothing with over flats SW 5.0 mm 32 mm length of thread, drill bit 5.0 mm

127-080-40	40 mm	127-080-90	90 mm
127-080-45	45 mm	127-080-91	95 mm
127-080-50	50 mm	127-080-92	100 mm
127-080-55	55 mm	127-080-93	105 mm
127-080-60	60 mm	127-080-94	110 mm
127-080-65	65 mm	127-080-95	115 mm
127-080-70	70 mm	127-080-96	120 mm
127-080-75	75 mm		
127-080-80	80 mm		
127-080-85	85 mm		




127-085-13	13 mm \varnothing	Washers for Cannulated Screws 3.5 mm
127-085-16	16 mm \varnothing	Washers for Cannulated Screws 4.5 mm
127-085-19	19 mm \varnothing	Washers for Cannulated Screws 7.0 mm




- 127-125-10 Guide wire 1.0 mm for Cannulated Screws 3.5 mm \varnothing 230 mm
- 127-125-16 Guide wire 1.6 mm for Cannulated Screws 4.5 mm \varnothing 230 mm
- 127-125-20 Guide wire 2.0 mm for Cannulated Screws 7.0 mm \varnothing 230 mm




- 127-160-27 Cannulated drill bit 2.7 mm for 3.5 mm Screws
- 127-160-35 Cannulated drill bit 3.5 mm for 3.5 mm Screws
- 127-160-32 Cannulated drill bit 3.2 mm for 4.5 mm Screws
- 127-160-45 Cannulated drill bit 4.5 mm for 4.5 mm Screws
- 127-160-50 Cannulated drill bit 5.0 mm for 7.0 mm Screws



- 127-168-35 Cannulated Tap for 3.5 mm Screws
- 127-168-45 Cannulated Tap for 4.5 mm Screws
- 127-168-70 Cannulated Tap for 7.0 mm Screws



- 127-173-35 Cannulated Screw Driver for 3.5 mm screws (fiber handle)
- 127-173-45 Cannulated Screw Driver for 4.5 mm screws (fiber handle)
- 127-173-70 Cannulated Screw Driver for 7.0 mm screws (fiber handle)



- 127-175-35 Cannulated Screw driver for 3.5 mm Screws with (T-handle)
- 127-175-45 Cannulated Screw driver for 4.5 mm Screws with (T-handle)
- 127-175-70 Cannulated Screw driver for 7.0 mm Screws with (T-handle)




127-200-35 Tap Sleeve 3.5 mm for 3.5 mm Screws
Drill Sleeve 2.5 mm



127-212-35 Tap Sleeve 2.7 mm for 3.5 mm Screws
Drill Sleeve 1.0 mm



100-115-08 Tap Sleeve 4.5 mm for 4.5 mm Screws
Drill Sleeve 3.2 mm



127-227-70 Tissue Protector Sleeve
for 7.0 mm Screws



127-245-20 Parallel Target Device for Guide Wires \varnothing 2.0 mm



Stainless Steel



- 128-002-16** Barbed Staples 8.7 x 16 x 2 mm
- 128-002-19** Barbed Staples 8.7 x 19 x 2 mm
- 128-002-22** Barbed Staples 8.7 x 22 x 2 mm
- 128-002-25** Barbed Staples 8.7 x 25 x 2 mm




	Diameter	Width	Length
128-015-22	1.5 mm	22 mm	19 mm
128-015-16	1.5 mm	16 mm	19 mm
128-025-22	2.5 mm	22 mm	19 mm
128-025-16	2.5 mm	16 mm	19 mm



128-033-48
24 x 33 x 4.8 mm

128-033-95
24 x 33 x 9.5 mm

128-033-14
24 x 33 x 14 mm



128-130-25 Staples Extractor




128-132-16 Staples Driver




128-134-19 Staples Inserter




- 128-250-06 6 mm
- 128-250-08 8 mm
- 128-250-11 11 mm
- 128-250-14 14 mm




128-270-18 Impactor




- 128-322-09 9.5 x 22 mm
- 128-322-16 16 x 22 mm
- 128-322-22 22 x 22 mm
- 128-322-28 28 x 22 mm




128-271-02 Impactor for 132-322-09 To 132-322-28
 128-271-04 Impactor for 132-340-16 To 132-355-15



Cat. No.:	Width	Length
128-340-16	16 mm	22 mm
128-340-22	22 mm	22 mm
128-340-28	28 mm	22 mm



Cat. No.:	Width A	Length B	Height C
128-348-05	24 mm	30 mm	5 mm
128-348-10	24 mm	30 mm	10 mm
128-348-15	24 mm	30 mm	15 mm



Cat. No.:	Width	Length	Height
128-355-05	24 mm	30 mm	5 mm
128-355-10	24 mm	40 mm	10 mm
128-355-15	24 mm	40 mm	15 mm