

Introduction


Indication/contraindication



PINE III Percutaneous Pedicle Screw System is designed for posterior thoracolumbar percutaneous fixation. The special "Self-tubular" design offers no extension loading or unloading and streamlines surgical technique. And the extreme long tab makes this system the capacity for obese patient.

- "Self-tubular" design

- Percutaneous fixation, minimize invasive to patient
- Ease of use, no extension loading or unloading



-360FAS pedicle screw design

- Mono-axial screw, screw head rotation along the vertical axis is allowed



-Dural-quadruple lead thread design

- Faster screw insertion
- Concordant screw insertion speed
- Respect pedicle anatomy, better pedicle purchase

