

A blue and white braided suture-tendon with a curved needle, shown against a light background. The suture is coiled in the center, with the needle pointing upwards and to the left. The text 'SUTUR'LINK®' is overlaid on the left side of the suture.

SUTUR'LINK®


Suture -Tendons
and Ligaments Reintegration

The logo for Teknimed, featuring a stylized blue arrow pointing right, followed by the word 'Teknimed' in a bold, italicized blue font.

Teknimed

Storage

SUTUR'LINK® is supplied sterile, ready to use in operating rooms. It must be kept in its original unopened packaging, in a dry, clean place away from light and at room temperature.



Product

SUTUR'LINK® is composed of a thread crimped on stainless steel needles. It is indicated to repair or reinforce ligaments, closure and/or ligation of soft tissues, and tuberosity's reinsertions.

- Thread: USP $\frac{3}{4}$, L = 900 mm
- Material: 87,5 % polyethylene / 12,5 % polypropylene
- Round needle: $\frac{3}{8}$ length 25 mm
- Triangular needle: $\frac{1}{2}$ length 40 mm

SUTUR'LINK® is sterilized with ethylene oxide. Single use. Do not re-sterilize. For any further information, please refer to the IFU.



DESIGNATION

Suture - Tendons and Ligaments Reintegration SUTUR'LINK®

REFERENCE

T362034

Distributed by

Manufactured by



Classe IIb



Head offices:

8, rue du Corps Franc-Pommiès
65500 VIC en BIGORRE (France)
Tél (33) 5 62 96 88 38
Fax (33) 5 62 96 28 72

Administration and supply-chain offices:

Z.I de la Herray
65500 VIC en BIGORRE (France)
Tél (33) 5 62 96 88 38
Fax (33) 5 62 96 28 72

Sales and engineering departments:

Z.I. de Montredon - 12, rue Apollo
31240 L'UNION (France)
Tél (33) 5 34 25 10 60
Fax (33) 5 34 25 27 39

